COVID-19 Antigen Rapid Test (Colloidal Gold) anterior nasal - self testing device
Instructions For Use
[PRODUCT NAME]

COVID-19 Antigen Rapid Test (Colloidal Gold) anterior nasal - self testing device

[SPECIFICATIONS AND COMPONENTS]

1 test/kit, 5 tests/kit, 10 tests/kit, 25 tests/kit.
 Test cassette, Extraction buffer, Sterile swab, Instructions for use

[INTENDED USE]

This product is used for the in vitro qualitative diagnosis of SARS-CoV-2 antigen in human nasal swab specimens. It is intended for personal use by persons as a rapid test method for coronavirus infection. The kit test results are for informational purposes only and are not a substitute for tests performed by healthcare professionals. It is suitable for users 0 years of age and older. Users under 14 years of age should be tested under supervision or with the assistance of an adult. Both symptomatic and asymptomatic individuals can be tested.

[TEST PRINCIPLE]

The COVID-19 Antigen Test (Colloidal Gold) is a qualitative membrane strip based on immunoassay for the detection of the nucleocapsid protein antigen from SARS-CoV-2 in swab specimen. If the sample contains the coronavirus antigen, both the test line (T) and control line (C) will appear and the result will be positive. If the sample does not contain coronavirus antigen or no coronavirus antigen is detected, the test line (T) will not appear, only the control line (C) appears and the result will be negative.

[STORAGE AND SHELF LIFE]


The test is stable for 2 years if all parts are stored in sealed pouches at +2°C to 50°C. The test should be stored in the sealed pouch until use.

[SAMPLE REQUIREMENTS]

Read the entire leaflet before performing the test.

Collection of samples:


- Please disinfect your hands thoroughly before use!
- Then open the package containing the sterile swab and take it out.
- Gently insert the entire soft end of the swab about 2 cm into one nostril.

- Using medium pressure, slowly rub the swab in a circular motion around the inside of your nostril 4 times for a total of 15 seconds.
- In the solution, rotate the swab and squeeze the side of the tube, repeating this three times.
- Repeat the same procedure with the same swab in the other nostril.

Sample handling: Mixing the sample with the extraction liquid

- Open the tube containing the extraction buffer.
- Dip the soft end of the sampling swab below the liquid level of the sample extraction buffer after sampling.
- In the solution, rotate the swab and squeeze the side of the tube, repeating this three times. The soft end of the sampling swab should be kept in the solution for at least 15 seconds.
- Then press the soft end of the swab three times against the inner wall of the tube.
- Dispose the swab in the waste bag.
- Close the tube containing the buffer solution with the dropper cap.

**Storage of samples:**


Drop the sample onto the test cassette as soon as possible, but no later than one hour after the previous process. If you do not perform the test for more than an hour after sampling, the test result will be considered invalid.

[TEST PROCEDURE]

- The test should be performed at room temperature (15°C - 30°C). Before performing the test, make sure that both the test cartridge and the extraction buffer are at room temperature.
 - Remove the test cassette from the sealed pouch and place the test cassette on a flat surface.
 - Place 2 drops of the sample-containing liquid prepared as described in the previous step into the sample hole (marked with the letter S).
 - Please set a timer to 15 minutes. The result can be read 15 minutes after dropping the sample-containing liquid on the cassette.
- The check line marked C may appear earlier. However, be sure to wait 15 minutes before reading. The results you read before this time are invalid. If no result is available after 20 minutes, it is recommended to repeat the test with a new test cassette.

**[INTERPRETATION OF RESULTS]**

- Positive (+):** Two lines appear (at C and T). One line should always appear in the control line area (C) and another line should appear in the test line area (T) (see Result 1).
- Negative (-):** A colored line appears at the control area (C). No colored line appears in the test line area (T) (see Result 2).
- Invalid:** The colored line does not appear at the control area (C) (see result 3). Check the instructions for use and repeat the test with a new kit. If the problem persists, stop using the test kit immediately and contact your local distributor.

A negative result only gives information about your current condition and does not mean that you do not carry the SARS-CoV-2 virus. The result obtained is for information only, SARS-CoV-2 infection can only be completely ruled out by a nucleic acid detection test (eg PCR test).

[DISPOSAL OF THE TEST KIT]

Place the test cassette, sample extraction buffer and used swab in a waste bag and seal the bag. Throw the bag in the bin and re-apply the hand sanitizer.

[RESTRICTIONS]

- The results of this test are for information only. The infection should be confirmed by a healthcare professional by comparing this test with other laboratory results, clinical symptoms, epidemiological data, and additional clinical data.
- In the early stages of infection, low virus concentrations may result in a negative test result.
- Test results depend on the quality of sample collection, processing, transport and storage. Errors can lead to false negative results. If contamination is not avoided during sample processing, a false positive result may be obtained.

[SENSITIVITY AND SPECIFICITY]

Method	RT-PCR		Total Results
	Positive	Negative	
COVID-19 Antigen Rapid Test	Positive	4	102
	Negative	496	500
Total Result	102	500	602

Diagnostic sensitivity: 96.1% (90.26%-98.92%)*
 Diagnostic specificity: 99.2% (97.86%-99.92%)*
 Overall agreement: 98.7% (97.40%-99.42%)*
 * 95% confidence interval.

[CROSS-REACTIONS]

No cross-reactivity was observed including HCoV-229E, HCoV-OC43, HCoV-NL63, MERS-CoV, HCoV-HKU1, Respiratory syncytial virus, Human enterovirus, Human rhinovirus, Human metapneumovirus, Mycoplasma pneumoniae, Parainfluenza, Adenovirus, Influenza B virus (Victoria line), H1N1 (2009) influenza virus, Influenza A virus H3N2, Avian influenza virus H7N9, Influenza B virus (Yamagata series), Neisseria meningitidis, Seasonal influenza A H1N1, Streptococcus pneumoniae, Staphylococcus aureus

[INTERFERING SUBSTANCES]

The test cassettes were tested for possible cross-reactions with the following chemicals and drugs. No cross-reactivity was observed including Aspirin, Ascorbic Acid, Ibuprofen, Bilirubin and Chloram-phenicol.

[PRECAUTIONS]

- This kit is intended for in vitro diagnostic use only. Please read these instructions carefully before testing.
- Use only the sample swab and extraction buffer provided in the package and do not substitute components from other kits for the specimen extracted provided in this kit.
- Strictly follow the instructions when performing the test.
- Positive and negative predictive values are highly dependent on prevalence. When disease prevalence is low and SARS-CoV-2 is not or only minimally active, positive test results are more likely to represent a false positive result. When disease prevalence is high, false negative results are more likely.
- Compared to a SARS-CoV-2 RT-PCR test, this test is less sensitive when used to detect patient samples within the first five days of symptom onset.
- Waste and excess samples generated during the test may be disposed of with household waste.
- The extraction buffer is not suitable for consumption.

[EFFECTIVE DATE AND VERSION]

Effective Date: 30.12.2021
 Version: 0

Note: Please refer to the table below to identify various symbols.

	Read the instruction Manual		Self-test
	Shelf life until		Tests per kit
	Lot		In vitro diagnostics
	Catalogue number		Attention
	Manufacturer		Date of Production
	Representative in the EU		Temperature limit
	Do not reuse		CE/ Directive 98/79/EC

Joinstar Biomedical Technology Co., Ltd.
 Address: 10th Floor, Administration Building, NO.519, XingGuo RD., Yuhang Economic and Technological Development Zone Hangzhou, Zhejiang, China, 311188
 E-mail: market@joinstar.cn
 Tel: 0086-571-89028135
 Fax: 0086-571-89028135

Lotus NL B.V.
 Address: Koninigin Julianaplein 10, 1e Verd, 2595AA, The Hague, Netherlands.
 E-mail: peter@lotsus.nl
 Tel: +31644168999

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