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Instructions for use:

1. Inspect the device, particularly the cap seal, for visible damage or signs of deterioration. If in doubt, replace the device.

2. Ensure that the device is dry and contains **no visible wetness or moisture**, and that the device is not warm to the touch. If in doubt, leave the device to cool with the cap removed for 15 minutes after a prior use.

3. Fold the indicator strip in half so that the yellow indicator faces inwards, then place it into the capsule of the device. The open end of the folded indicator should face the open end of the capsule.

4. Check that the indicator does not protrude from the capsule.

5. Close the device, ensuring that the cap is sealed. The device has two lugs to serve as a 'stop'.

6. Seal the device into a paper / plastic sterilization pouch of dimensions (approx) 110mm x 230mm or 90mm x 230mm, or 133mm x 254mm and place it into the centre of the sterilizer.

7. Process the device in an otherwise empty chamber using the wrapped instrument cycle at 132°C - 137°C for up to 3.5 minutes.

8. At the end of the cycle, remove the device from the sterilizer within 10 minutes.

Note:

Caution is required as the device may be hot.

BROWNE TST Helix Steam Penetration Test Device

This device satisfies the requirement to perform a daily steam penetration test given in:

- **EN ISO 17665-1:2006** Specifies requirements for the development, validation and routine control of a moist heat sterilization process for medical devices

- **DB 2002(06)** Device Bulletin, Benchmark Steam Sterilizers – Guidance on purchase, Operation and maintenance, Medical Devices Agency, clause 5.3.1

- **SN 2002(24)** Safety Notice, Steam penetration tests in vacuum benchtop sterilizers, medical Devices Agency

- **DB 9804** Device bulletin, The validation and periodic testing of benchtop vacuum steam sterilizers, medical Devices Agency, clause 5.1


Note:

Some of the above documents are available on the UK Medical Devices Agency's website, www.medical-devices.gov.uk


Description:

The device consists of a plastic tubular device with an indicator receptacle at the end and is supplied with indicator strips. After 250 indicators have been used, the device must be replaced.


9. Remove the indicator from the capsule of the device and check the indicator result:
Yellow = Unprocessed

Presence of Yellow/Brown/Green = Fail

Blue/Purple = Pass

10. Record the indicator result in the sterilizer log book, along with the cycle details, and retain the indicator if required. Each indicator is provided with a self adhesive backing to facilitate attachment to sterilizer records.

11. If the result of the indicator is a PASS result, i.e. blue indicator result, the air removal stage of the sterilizer is working properly, and wrapped medical devices may be processed in the sterilizer (see instructions for use for the sterilizer).

12. Store the device with the cap removed to facilitate drying of the tubing.

IMPORTANT:

The indicators must be stored in a cool, dark and dry place before and after use.

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2367-20r02
April 2008

Use:

In order to sterilize products and devices by steam, steam must be able to penetrate through and contact all surfaces. Items such as packaging, pouches, narrow or hollow instruments are not easy to penetrate with steam. For this reason, steam sterilizers are fitted with vacuum pumps to remove the air and thus allow steam to penetrate into the load. This device is designed to test the steam penetration ability of the sterilizer, by using a device known to represent a defined challenge to the sterilizer.

Frequency of use:

To comply with the requirements of the regulatory documents listed above, a steam penetration test must be conducted at the commencement of each day that the sterilizer is to be used. This 'verifies' the suitability of the sterilizer's air removal stage for the day ahead.

IMPORTANT:

- The device has been developed and validated solely for use in Type B Steam Sterilizers
- Use of this device in any other sterilizer, or use of this device with any other type of indicator, may give dangerously misleading results
- Failure to replace the device after 250 uses could lead to failure of the device and dangerously misleading results

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